



Model MI-30

30 Watt Convection Cooled Illuminator

Standard Features:

- 10,000 Footcandles
- Intense, Cold Illumination
- Convection-Cooled
- Sleek, Stackable Housing
- 3-Year Warranty
- UL, C-UL

Accessories & Options:

- Optional Single Gooseneck Light Guide
- Optional Filter Holder & IR Filter





Technical Information

Input:	115 V AC
Lamp:	30 Watt, 10.8 V, Type EKZ qtz halogen
Color Temp.:	3100° Kelvin
Lamp Life:	200 hrs. full intensity
Lamp Output:	10,000 Foot Candles
Fiber Optic Interface:	0.59"
Noise Level:	Silent, No Fan, No Vibration
Safety Certifications:	UL, C-UL
Dimensions:	6.0" x 8.5" x 5.3"; 7.25 lbs.

Options

MI-30:	30 Watt, 115 V AC Illuminator
MI-30F:	30 Watt, 115 V AC Illuminator w/ Filter Holder & IR Filter
MI-30SG:	30 Watt, 115 V AC Illuminator w/ BGG1823M Single Gooseneck, (1) LH759 Focusing Lense
MI-30SGF:	30 Watt, 115 V AC Illuminator w/ BGG1823M Single Gooseneck, (1) LH759 Focusing Lense, Filter Holder & IR Filter

MI-30 Dimensions:



[mm]
in

Dolan-Jenner
industries

Dolan-Jenner Industries
159 Swanson Rd.
Boxborough, MA 01719 USA

Phone: 800.833.4237
Fax: 978.264.0292
Email: support@dolan-jenner.com

1.800.83.FIBER

WWW.DOLAN-JENNER.COM