



A645NK41



700009570036



700010730039


A645NK41 Series - Nikon Microscope Ringlight A

Mounting: Threaded Locking Ring (Included) Working Distance: 1.25in. – 5.0in.

Model	Part #	Cable Length (in.)	End Tip Dia. (in.)	Light Source Adapter .590 (15mm) Input	Light Source Adapter 1.00 (25mm) Input	Application
A645NK41	700008560044	41	.588	Direct Fit	SX-13	SMZ645, SMZ660, SMZ660V, SMZ1000

700009570036 Series - Nikon Microscope Ringlight C

Mounting: Threaded Locking Ring (Included) Working Distance: 1.25in. – 5.0in.


Model	Part #	Cable Length (in.)	End Tip Dia. (in.)	Light Source Adapter .590 (15mm) Input	Light Source Adapter 1.00 (25mm) Input	Optional Polarizer
700009570036	700009570036	36	.588	Direct Fit	SX-13	LP-PLZ 700009730000 

70001073 Series - Nikon Microscope Large Diameter Ring Light


Mounting: Threaded Locking Ring (Included) Working Distance: 2.0in. – 6.0in.

Model	Part #	Cable Length (in.)	End Tip Dia. (in.)	Light Source Adapter .590 (15mm) Input	Light Source Adapter 1.00 (25mm) Input	Application
700010730039	700010730039	39	.588	Direct Fit	SX-13	SMZ800(w/ .5 X 1.5 and 2 objective lenses), SMZ1000(w/ .5 X 1.5 and 2 objective lenses), SMZU(w/.5 X 1.5 objective lens), SMZ10A(w/.5 X 1.5 objective lens), SMZ1500(w/.5 X 1 objective lens)
700010730060	700010730060	60	.588	Direct Fit	SX-13	SMZ800(w/ .5 X 1.5 and 2 objective lenses), SMZ1000(w/ .5 X 1.5 and 2 objective lenses), SMZU(w/.5 X 1.5 objective lens), SMZ10A(w/.5 X 1.5 objective lens), SMZ1500(w/.5 X 1 objective lens)


A645NK41 Dimensions:



700009570036 Dimensions:



70001073 Dimensions:



Note: Ring lights available in variety of lengths. 36 inch version shown in diagram above. Please contact Dolan-Jenner for availability on length variations. Please note two standard end ferrules available along with many custom designs.

Dolan-Jenner
industries

Dolan-Jenner Industries
159 Swanson Rd.
Boxborough, MA 01719 USA

Phone: 800.833.4237
Fax: 978.264.0292
Email: support@dolan-jenner.com

1.800.83.FIBER

WWW.DOLAN-JENNER.COM